

anti- RIN3 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab07307 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Ras effector protein that functions as a guanine nucleotide exchange(GEF) for RAB5B and RAB31, by exchanging bound GDP for free GTP. Required for normal RAB31 function.

Immunogen information

| | |
|--------------|--------------------------------|
| Immunogen: | Ras and Rab interactor 3 |
| Synonyms: | Ras and Rab interactor 3, RIN3 |
| Observed MW: | 108 kDa |
| UniprotID : | Q8TB24 |

Application


| | |
|--------------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution:WB: | 1:200-1:1000 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07307(RIN3 antibody) at dilution of 1:300