

anti- RND3 antibody

Product Information

Catalog No.:	FNab07334
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2F0
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Binds GTP but lacks intrinsic GTPase activity and is resistant to Rho-specific GTPase-activating proteins.

Immunogen information

Immunogen:	Rho family GTPase 3
Synonyms:	ARHE, memB, Protein MemB, Rho family GTPase 3, Rho8, RhoE, RND3
Observed MW:	29 kDa
Uniprot ID :	P61587

Application

Reactivity:	Human, Pig
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab07334(RND3 Antibody) at dilution of 1:1000