

anti- SCML2 antibody

Product Information

Catalog No.:	FNab07640
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Putative Polycomb group(PcG) protein. PcG proteins act by forming multiprotein complexes, which are required to maintain the transcriptionally repressive state of homeotic genes throughout development(By similarity).

Immunogen information

Immunogen:	sex comb on midleg-like 2(Drosophila)
Synonyms:	None
Observed MW:	70 kDa
Uniprot ID :	Q9UQR0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab07640(SCML2 Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab07640(SCML2 Antibody) at dilution of 1:300