

# anti- SDHB antibody

## **Product Information**

Catalog No.:	FNab07670
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

## Background

Iron-sulfur protein(IP) subunit of succinate dehydrogenase(SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone(coenzyme Q).

### **Immunogen information**

Immunogen:	succinate dehydrogenase complex, subunit B, iron sulfur(Ip)
Synonyms:	IP, PGL4, SDH, SDH1, SDHB, SDHIP
Observed MW:	32 kDa
UniprotID :	P21912


## Application

Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, WB, IHC, IP	
Recommended dilution: WB: 1:500-1:2000; IHC: 1:100-1:500		
Image:		

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab07670(SDHB Antibody) at dilution of 1:200

IP Result of anti-SDHB (IP:FNab07670, 4ug; Detection:FNab07670 1:500) with HEK-293 cells lysate 2800ug.




#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>