

anti- SEC61B-Specific antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab07691 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

None

Immunogen information

| | |
|--------------|--------------------|
| Immunogen: | Sec61 beta subunit |
| Synonyms: | None |
| Observed MW: | 10 kDa |
| UniprotID : | P60468 |

Application


| | |
|-----------------------|--|
| Reactivity: | Human, Mouse |
| Tested Application: | ELISA, WB, IHC, IF |
| Recommended dilution: | WB: 1:1000-1:4000; IF: 1:20-1:200; IHC: 1:20-1:200 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using FNab07691(SEC61B-Specific Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab07691(SEC61B-Specific antibody) at dilution of 1:500