

anti- SEPT4 antibody

Product Information

Catalog No.:	FNab07723
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Filament-forming cytoskeletal GTPase(By similarity). May play a role in cytokinesis(Potential). May play a role in platelet secretion. Isoform ARTS, but not the other isoforms, is required for the induction of cell death mediated by TGF-beta and by other apoptotic stimuli.

Immunogen information

Immunogen:	septin 4
Synonyms:	ARTS, PNUTL2, SEP4
Observed MW:	30 kDa
Uniprot ID :	O43236

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



A375 cells were subjected to SDS PAGE followed by western blot with FNab07723(SEPT4 antibody) at dilution of 1:300