

anti- SEPT8 antibody

Product Information

Catalog No.:	FNab07728
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9H10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Filament-forming cytoskeletal GTPase(By similarity). May play a role in cytokinesis(Potential). May play a role in platelet secretion.

Immunogen information

Immunogen:	septin 8
Synonyms:	KIAA0202
Observed MW:	53 kDa
Uniprot ID :	Q92599

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07728(SEPT8 Antibody) at dilution of 1:1000