

anti- SFTPC antibody

Product Information

Catalog No.:	FNab07798
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Pulmonary surfactant associated proteins promote alveolar stability by lowering the surface tension at the air-liquid interface in the peripheral air spaces.

Immunogen information

Immunogen:	surfactant protein C
Synonyms:	PSP C, SFTP2, SFTPC, SMDP2, SP C, SP5, surfactant protein C
Observed MW:	21 kDa
Uniprot ID :	P11686

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human lung cancer using FNab07798(SFTPC antibody) at dilution of 1:50



A549 cells were subjected to SDS PAGE followed by western blot with FNab07798(SFTPC antibody) at dilution of 1:1000