

## anti- SH3BP1 antibody

### Product Information

Catalog No.:	FNab07829
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Binds differentially to the SH3 domains of certain proteins of signal transduction pathways. This protein binds preferentially to the ABL1 proto-oncogene, SRC and GRB2. Shows GAP activity for Rac-related proteins but not for Rho-or Ras-related proteins. It inhibits PDGF-induced membrane ruffling mediated by Rac(By similarity).

### Immunogen information

Immunogen:	SH3-domain binding protein 1
Synonyms:	None
Observed MW:	76-80 kDa,51kDa
UniprotID :	Q9Y3L3

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:WB:	1:500-1:2000; IF: 1:20-1:200


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Image:



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab07829(SH3BP1 antibody) at dilution of 1:500

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)