

## anti- SKIV2L2 antibody

### Product Information

Catalog No.:	FNab07896
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

### Immunogen information

Immunogen:	superkiller viralicidic activity 2-like 2 (S. cerevisiae)
Synonyms:	DOB1, KIAA0052, Mtr4
Observed MW:	123 kDa
Uniprot ID :	P42285

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human lung cancer using FNaB07896( SKIV2L2 Antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNaB07896( SKIV2L2 Antibody) at dilution of 1:1000