

## anti- SLC22A18AS antibody

### Product Information

Catalog No.:	FNab07922
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	solute carrier family 22(organic cation transporter), member 18 antisense
Synonyms:	BWR1B, BWSCR1B, ORCTL2S, SLC22A1LS
Observed MW:	None
UniprotID :	Q8N1D0

### Application


Reactivity:	Human
Tested Application:	ELISA, IP
Recommended dilution:IP:	1:200-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-SLC22A18AS (IP:FNab07922, 3ug; Detection:FNab07922 1:300) with HepG2 cells lysate 3600ug.