

## anti- SLC36A4 antibody

### Product Information

Catalog No.:	FNab07955
Size:	100µ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Functions as a sodium-independent electroneutral transporter for tryptophan, proline and alanine. Inhibited by sarcosine.

### Immunogen information

Immunogen:	solute carrier family 36(proton/amino acid symporter), member 4
Synonyms:	Synonyms PAT4, SLC36A4
Calculated MW:	41-66 kDa
Uniprot ID :	Q6YBV0

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



A549 cells were subjected to SDS PAGE followed by western blot with FNab07955(SLC36A4 antibody) at dilution of 1:1000