

## anti- SLITRK6 antibody

### Product Information

Catalog No.:	FNab07984
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

### Background

Regulator of neurite outgrowth required for normal hearing and vision.

### Immunogen information

Immunogen:	SLIT and NTRK-like family, member 6
Synonyms:	SLIT and NTRK like protein 6, SLITRK6
Observed MW:	95-100 kDa
Uniprot ID :	Q9H5Y7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IF: 1:10-1:100
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunofluorescent analysis of HepG2 cells, using SLITRK6 antibody FNab07984 at 1:25 dilution. Blue pseudocolor = DAPI (fluorescent DNA dye).



A549 cells were subjected to SDS PAGE followed by western blot with FNab07984(SLITRK6 antibody) at dilution of 1:500