

---

## anti- Somatostatin antibody

### Product Information

Catalog No.:	FNab08109
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	somatostatin
Synonyms:	SMST, somatostatin, SST
Observed MW:	17 kDa
Uniprot ID :	P61278

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:500
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab08109(SST Antibody) at dilution of 1:200



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab08109(SST antibody) at dilution of 1:1000