

anti- SPARCL1 antibody

Product Information

Catalog No.:	FNab08154
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	SPARC-like 1(hevin)
Synonyms:	None
Observed MW:	130-140 kDa
UniprotID :	Q14515

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using FNab08154 (SPARCL1 Antibody) at dilution of 1:50



IP Result of anti-SPARCL1 (IP:FNab08154, 4ug; Detection:FNab08154 1:300) with Raji cells lysate 1440ug.



human heart tissue were subjected to SDS PAGE followed by western blot with FNab08154(SPARCL1 antibody) at dilution of 1:400