

anti- ST8SIA1 antibody

Product Information

Catalog No.:	FNab08275
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Involved in the production of gangliosides GD3 and GT3 from GM3; gangliosides are a subfamily of complex glycosphingolipids that contain one or more residues of sialic acid.

Immunogen information

Immunogen:	ST8 alpha-N-acetyl-neuraminide alpha-2,8-sialyltransferase 1
Synonyms:	SIAT8, SIAT8A
Observed MW:	43 kDa
UniprotID :	Q92185

Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab08275(ST8SIA1 Antibody) at dilution of 1:1000