

anti- STRA6 antibody

Product Information

Catalog No.:	FNab08349
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May act as a high-affinity cell-surface receptor for the complex retinol-retinol binding protein(RBP/RBP4). Acts by removing retinol from RBP/RBP4 and transports it across the plasma membrane, where it can be metabolized. This mechanism does not depend on endocytosis. Binds to RBP/RBP4 with high affinity. Increases cellular retinol uptake from the retinol-RBP complex(By similarity).

Immunogen information

Immunogen:	stimulated by retinoic acid gene 6 homolog(Mouse)
Synonyms:	None
Observed MW:	None
Uniprot ID :	Q9BX79

Application

Reactivity:	Human
Tested Application:	ELISA
Recommended dilution:	None

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com