

anti- TAB1 antibody

Product Information

Catalog No.:	FNab08464
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be an important signaling intermediate between TGFβ receptors and MAP3K7/TAK1. May play an important role in mammalian embryogenesis.

Immunogen information

Immunogen:	mitogen-activated protein kinase kinase kinase 7 interacting protein 1
Synonyms:	MAP3K7IP1
Observed MW:	55 kDa
Uniprot ID :	Q15750

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889


www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse lung tissue slide using FNab08464(TAB1 antibody) at dilution of 1:200



Immunofluorescent analysis of U251MG cells, using TAB1 antibody FNab08464 at 1:50 dilution.



K-562 cells were subjected to SDS PAGE followed by western blot with FNab08464(TAB1 antibody) at dilution of 1:600