

## anti- TCTN2 antibody

### Product Information

Catalog No.:	FNab08568
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Component of the tectonic-like complex, a complex localized at the transition zone of primary cilia and acting as a barrier that prevents diffusion of transmembrane proteins between the cilia and plasma membranes. Required for hedgehog signaling transduction(By similarity).

### Immunogen information

Immunogen:	tectonic family member 2
Synonyms:	C12orf38, TECT2
Observed MW:	77 kDa
Uniprot ID :	Q96GX1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IP, IF, IHC
Recommended dilution:	IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

1


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemical of paraffin-embedded human liver using FNab08568(TCTN2 antibody) at dilution of 1:100



IP Result of anti-TCTN2 (IP: FNab08568, 4ug; Detection: FNab08568 1:500) with mouse lung tissue lysate 4000ug.