

anti- TCTN3 antibody

Product Information

Catalog No.:	FNab08569
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Part of the tectonic-like complex which is required for tissue-specific ciliogenesis and may regulate ciliary membrane composition(By similarity). May be involved in apoptosis regulation. Necessary for signal transduction through the sonic hedgehog(Shh) signaling pathway.

Immunogen information

Immunogen:	tectonic family member 3
Synonyms:	C10orf61, TCTN3, TECT3, Tectonic 3, tectonic family member 3, UNQ1881/PRO4324
Observed MW:	66 kDa
UniprotID :	Q6NUS6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human lung using FNab08569(TCTN3 antibody) at dilution of 1:50



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08569(TCTN3 Antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com