

anti- TFF2 antibody

Product Information

Catalog No.:	FNab08629
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Inhibits gastrointestinal motility and gastric acid secretion. Could function as a structural component of gastric mucus, possibly by stabilizing glycoproteins in the mucus gel through interactions with carbohydrate side chains(By similarity).

Immunogen information

Immunogen:	trefoil factor 2
Synonyms:	SML1, SP, Spasmolysin, Spasmolytic polypeptide, TFF2, trefoil factor 2
Observed MW:	18kd
Uniprot ID :	Q03403

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:1000; IP: 1:500-1:2000; IHC: 1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemical of paraffin--embedded mouse liver using FNab08629(TFF2 antibody) at dilution of 1:100



human stomach tissue were subjected to SDS PAGE followed by western blot with FNab08629(TFF2 antibody) at dilution of 1:300