

anti- TFPT antibody

Product Information

Catalog No.:	FNab08635
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Appears to promote apoptosis in a p53/TP53-independent manner. Putative regulatory component of the chromatin remodeling INO80 complex which is involved in transcriptional regulation, DNA replication and probably DNA repair.

Immunogen information

Immunogen:	TCF3(E2A) fusion partner(in childhood Leukemia)
Synonyms:	INO80F
Observed MW:	29 kDa
UniprotID :	P0C1Z6

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human tonsillitis slide using FNab08635(TFPT Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab08635(TFPT antibody) at dilution of 1:600