

anti- TMEM127 antibody

Product Information

Catalog No.:	FNab08750
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Controls cell proliferation acting as a negative regulator of TOR signaling pathway mediated by mTORC1. May act as a tumor suppressor.

Immunogen information

Immunogen:	transmembrane protein 127
Synonyms:	None
Observed MW:	None
UniprotID :	O75204

Application


Reactivity:	Human
Tested Application:	ELISA, IP, IF, IHC
Recommended dilution:	IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human heart slide using FNab08750(TMEM127 antibody) at dilution of 1:50



IP Result of anti-TMEM127 (IP:FNab08750, 3ug; Detection:FNab08750 1:500) with HeLa cells lysate 3000ug.