

## anti- TMEM213 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab08771   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

None

### Immunogen information

|              |                           |
|--------------|---------------------------|
| Immunogen:   | transmembrane protein 213 |
| Synonyms:    | None                      |
| Observed MW: | None                      |
| Uniprot ID : | A2RRL7                    |

### Application

|                            |            |
|----------------------------|------------|
| Reactivity:                | Human      |
| Tested Application:        | ELISA, IHC |
| Recommended dilution: IHC: | 1:20-1:200 |
| Image:                     |            |


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab08771(TM213 antibody) at dilution of 1:50