

## anti- TOP1 antibody

### Product Information

Catalog No.:	FNab08865
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

This gene encodes a DNA topoisomerase, an enzyme that controls and alters the topologic states of DNA during transcription. This enzyme catalyzes the transient breaking and rejoining of a single strand of DNA which allows the strands to pass through one another, thus altering the topology of DNA. This gene is localized to chromosome 20 and has pseudogenes which reside on chromosomes 1 and 22.

### Immunogen information

Immunogen:	topoisomerase (DNA) I
Synonyms:	DNA topoisomerase 1, DNA topoisomerase I, TOP1, TOPI, Top-I, topoisomerase (DNA) I
Observed MW:	100 kDa
UniprotID :	P11387

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IP

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IP: 1:50 - 1:100

Image:



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab08865(TOP1 antibody) at dilution of 1:50



L02 cells were subjected to SDS PAGE followed by western blot with FNab08865(TOP1 antibody) at dilution of 1:1000

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)