

## anti- TPPII antibody

### Product Information

Catalog No.:	FNab08896
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4A3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	tripeptidyl peptidase II
Synonyms:	FLJ40359, TPP 2, TPP II, TPP2, TPPII, Tripeptidyl aminopeptidase, Tripeptidyl peptidase 2, Tripeptidyl peptidase II
Observed MW:	135 kDa
Uniprot ID :	P29144

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human testis using FNab08896(TPP2 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08896(TPP2 antibody) at dilution of 1:1000