

anti- TRERF1 antibody

Product Information

Catalog No.:	FNab08962
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Binds DNA and activates transcription of CYP11A1. Interaction with CREBBP and EP300 results in a synergistic transcriptional activation of CYP11A1.

Immunogen information

Immunogen:	transcriptional regulating factor 1
Synonyms:	TRERF1, dJ139D8.5, FLJ21198, Zinc finger protein rapa, RP1 139D8.5
Observed MW:	132 kDa
UniprotID :	Q96PN7

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab08962(TRERF1 antibody) at dilution of 1:300