

anti- TRIM39 antibody

Product Information

Catalog No.:	FNab08984
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

E3 ubiquitin-protein ligase. May facilitate apoptosis by inhibiting APC/C-Cdh1-mediated poly-ubiquitination and subsequent proteasome-mediated degradation of the pro-apoptotic protein MOAP1.

Immunogen information

Immunogen:	tripartite motif-containing 39
Synonyms:	TRIM39, RNF23, TRIM39B, RING finger protein 23, TFP
Observed MW:	45 kDa, 56 kDa
Uniprot ID :	Q9HCM9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab08984(TRIM39 antibody) at dilution of 1:400