

## anti- TRIM44 antibody

### Product Information

Catalog No.:	FNab08986
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May play a role in the process of differentiation and maturation of neuronal cells(By similarity).  
May regulate the activity of TRIM17.

### Immunogen information

Immunogen:	tripartite motif-containing 44
Synonyms:	TRIM44, tripartite motif containing 44, Protein DIPB, MC7, HSA249128
Observed MW:	50 kDa
UniprotID :	Q96DX7

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human gliomas using FNab08986(TRIM44 antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08986(TRIM44 antibody) at dilution of 1:500