

anti- TRIM44 antibody

Product Information

Catalog No.:	FNab08987
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9A2
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May play a role in the process of differentiation and maturation of neuronal cells(By similarity).
May regulate the activity of TRIM17.

Immunogen information

Immunogen:	tripartite motif-containing 44
Synonyms:	TRIM44, tripartite motif containing 44, Protein DIPB, MC7, HSA249128
Observed MW:	50 kDa
Uniprot ID :	Q96DX7

Application

Reactivity:	Human
Tested Application:	ELISA, WB,IF
Recommended dilution:	WB: 1:500-1:5000; IF: 1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab08987(TRIM44 Antibody) at dilution of 1:1000