

## anti- TRMT6 antibody

### Product Information

Catalog No.:	FNab09008
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Substrate-binding subunit of tRNA(adenine-N(1)-)-methyltransferase, which catalyzes the formation of N(1)-methyladenine at position 58(m1A58) in initiator methionyl-tRNA.

### Immunogen information

Immunogen:	tRNA methyltransferase 6 homolog( <i>S. cerevisiae</i> )
Synonyms:	CGI 09, GCD10, Gcd10p, KIAA1153, TRM6, TRMT6, tRNA(m1A58)MTase subunit TRM6
Observed MW:	56 kDa
UniprotID :	Q9UJA5

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab09008 (TRMT6 Antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab09008 (TRMT6 antibody) at dilution of 1:500