

anti- TRPT1 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab09027 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Catalyzes the last step of tRNA splicing, the transfer of the splice junction 2'-phosphate from ligated tRNA to NAD to produce ADP-ribose 1"-2" cyclic phosphate.

Immunogen information

| | |
|--------------|---|
| Immunogen: | tRNA phosphotransferase 1 |
| Synonyms: | TRNA 2 phosphotransferase 1, tRNA phosphotransferase 1, TRPT1 |
| Observed MW: | 98 kDa |
| Uniprot ID : | Q86TN4 |

Application

| | |
|---------------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, IHC |
| Recommended dilution:IHC: | 1:20-1:200 |
| Image: | |


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab09027(TRPT1 Antibody) at dilution of 1:50