

---

## TSPAN2 antibody

### Product Information

Catalog No.:	FNab09059
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


May play a role in signalling in oligodendrocytes in the early stages of their terminal differentiation into myelin-forming glia and may also function in stabilizing the mature sheath.

### Immunogen information


Immunogen:	tetraspanin 2
Synonyms:	6330415F13Rik, FLJ12082, NET3, Tetraspan NET 3, tetraspanin 2, TSN2, TSPAN 2, TSPAN2
Observed MW:	30 kDa
Uniprot ID :	O60636

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-TSPAN2 (IP:FNab09059, 4ug;  
Detection:FNab09059 1:500) with HL-60 cells  
lysate 960ug.



HL-60 cells were subjected to SDS PAGE  
followed by western blot with  
FNab09059(TSPAN2 antibody) at dilution of  
1:500