

anti- TSPYL4 antibody

Product Information

Catalog No.:	FNab09069
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	TSPY-like 4
Synonyms:	DJ486I3.2, KIAA0721, TSPY like 4, TSPY like protein 4, TSPYL4
Observed MW:	45 kDa
Uniprot ID :	Q8VD63

Application

Reactivity:	Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09069(TSPYL4 Antibody) at dilution of 1:1000