

anti- TSTA3 antibody

Product Information

Catalog No.:	FNab09075
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Catalyzes the two-step NADP-dependent conversion of GDP-4-dehydro-6-deoxy-D-mannose to GDP-fucose, involving an epimerase and a reductase reaction.

Immunogen information

Immunogen:	tissue specific transplantation antigen P35B
Synonyms:	TSTA3, Protein FX, SDR4E1, FX, P35B, GDP-L-fucose synthase, GDP-4-keto-6-deoxy-D-mannose-3,5-epimerase-4-reductase, Red cell NADP(H)-binding protein, Short-chain dehydrogenase/reductase family 4E member 1
Observed MW:	40 kDa, 68 kDa
Uniprot ID :	Q13630

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000

1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

Image:



human colon tissue were subjected to SDS PAGE followed by western blot with FNab09075(TSTA3 antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com