

anti- TTC30A antibody

Product Information

Catalog No.:	FNab09079
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Required for polyglutamylation of axonemal tubulin. Plays a role in anterograde intraflagellar transport (IFT), the process by which cilia precursors are transported from the base of the cilium to the site of their incorporation at the tip.

Immunogen information

Immunogen:	tetratricopeptide repeat domain 30A
Synonyms:	FLJ13946, FLJ77601, IFT70, TPR repeat protein 30A, TTC30A
Observed MW:	76 kDa
Uniprot ID :	Q86WT1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000; IP: 1:200-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



rat testis tissue were subjected to SDS PAGE followed by western blot with FNab09079(TTC30A Antibody) at dilution of 1:1000