

anti- TUBGCP3 antibody

Product Information

Catalog No.:	FNab09102
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Immunogen information

Immunogen:	tubulin, gamma complex associated protein 3
Synonyms:	GCP 3, GCP3, h104p, hGCP3, hGrip104, hSpc98, SPBC98, Spc98p, TUBGCP3
Observed MW:	105 kDa
UniprotID :	Q96CW5

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver injury using FNab09102(TUBGCP3 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09102(TUBGCP3 antibody) at dilution of 1:1000