

## anti- TUBGCP4 antibody

### Product Information

Catalog No.:	FNab09103
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.

### Immunogen information

Immunogen:	tubulin, gamma complex associated protein 4
Synonyms:	76P, FLJ14797, GCP 4, GCP4, h76p, hGCP4, hGrip76, TUBGCP4
Observed MW:	76 kDa
Uniprot ID :	Q9UGJ1

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human placenta using FNab09103(TUBGCP4 antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab09103(TUBGCP4 antibody) at dilution of 1:500