

## anti- UBAP1 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab09149   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

Component of the ESCRT-I complex, a regulator of vesicular trafficking process. Binds to ubiquitinated cargo proteins and is required for the sorting of endocytic ubiquitinated cargos into multivesicular bodies(MVBs). Plays a role in the proteasomal degradation of ubiquitinated cell-surface proteins, such as EGFR and BST2.

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | ubiquitin associated protein 1                                  |
| Synonyms:    | NAG20, UAP, UBAP, UBAP 1, UBAP1, ubiquitin associated protein 1 |
| Observed MW: | 55 kDa  |
| Uniprot ID : | Q9NZ09  |

### Application

|                       |   |
|-----------------------|---|
| Reactivity:           | Human, Mouse, Rat   |
| Tested Application:   | ELISA, WB, IHC, IF, IP  |
| Recommended dilution: | WB: 1:500-1:2000; IF: 1:10-1:100; IP: 1:200-1:1000; IHC: 1:20-1:200 |

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)


Image:



Immunohistochemistry of paraffin-embedded human lung cancer using FNab09149(UBAP1 antibody) at dilution of 1:50



IP Result of anti-UBAP1 (IP:FNab09149, 4ug; Detection:FNab09149 1:500) with mouse brain tissue lysate 4000ug.



human heart tissue were subjected to SDS PAGE followed by western blot with FNab09149(UBAP1 antibody) at dilution of 1:500