

## anti- UBE2DNL antibody

### Product Information

Catalog No.:	FNab09170
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	ubiquitin-conjugating enzyme E2D N-terminal like pseudogene
Synonyms:	Putative ubiquitin-conjugating enzyme E2 D2-like protein, Ubiquitin carrier protein D2-like, Ubiquitin-conjugating enzyme E2D N-terminal-like, Ubiquitin-protein ligase D2-like, UBE2D2L
Observed MW:	23 kDa
Uniprot ID :	Q8IWF7

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution: IHC:	1:20-1:200
Image:	

1


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemical of paraffin-embedded human testis using FNab09170(UBE2DNL antibody) at dilution of 1:50