

anti- UBE2Q2 antibody

Product Information

Catalog No.:	FNab09182
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-48'-linked polyubiquitination.

Immunogen information

Immunogen:	ubiquitin-conjugating enzyme E2Q family member 2
Synonyms:	UBE2Q2, UBCi, Ubiquitin protein ligase Q2, DKFZp762C143, Ubiquitin carrier protein Q2
Observed MW:	43 kDa
Uniprot ID :	Q8WVN8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IP: 1:500-1:1000
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



IP result of anti-UBE2Q2(IP:FNab09182, 4ug;
Detection:FNab09182 1:300) with HEK-293T
cells lysate 1320 ug



A431 cells were subjected to SDS PAGE followed
by western blot with FNab09182(UBE2Q2
antibody) at dilution of 1:400