

## anti- UBE2Z antibody

### Product Information

Catalog No.:	FNab09188
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis.

### Immunogen information

Immunogen:	ubiquitin-conjugating enzyme E2Z
Synonyms:	UBE2Z, Ubiquitin carrier protein Z, FLJ13855, HOYS7, USE1
Observed MW:	25 kDa, 39 kDa
UniprotID :	Q9H832

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:100
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded mouse kidney using FNab09188(UBE2Z antibody) at dilution of 1:100



human testis tissue were subjected to SDS PAGE followed by western blot with FNab09188(UBE2Z antibody) at dilution of 1:1000

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)