

## anti- UCMA antibody

### Product Information

Catalog No.:	FNab09225
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May be involved in the negative control of osteogenic differentiation of osteochondrogenic precursor cells in peripheral zones of fetal cartilage and at the cartilage-bone interface.

### Immunogen information

Immunogen:	upper zone of growth plate and cartilage matrix associated
Synonyms:	C10orf49, Gla rich protein, GRP, UCMA, Ucma C
Observed MW:	33 kDa
Uniprot ID :	Q8WVF2

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09225(UCMA Antibody) at dilution of 1:1000