

## anti- UNKL antibody

### Product Information

Catalog No.:	FNab09266
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May participate in a protein complex showing an E3 ligase activity regulated by RAC1. Ubiquitination is directed towards itself and possibly other substrates, such as SMARCD2/BAF60b. Intrinsic E3 ligase activity has not been proven.

### Immunogen information

Immunogen:	unkempt homolog(Drosophila)-like
Synonyms:	C16orf28, FLJ12623, FLJ23360, KIAA0734, UNKL, ZC3H5L, ZC3HDC5L
Observed MW:	35 kDa, 45 kDa
Uniprot ID :	Q9H9P5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)

Image:



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09266(UNKL antibody) at dilution of 1:1500