

anti- UPB1 antibody

Product Information

Catalog No.:	FNab09269
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Converts N-carbamoyl-beta-aminoisobutyrate and N-carbamoyl-beta-alanine(3-ureidopropanoate) to, respectively, beta-aminoisobutyrate and beta-alanine, ammonia and carbon dioxide.

Immunogen information

Immunogen:	ureidopropionase, beta
Synonyms:	Beta alanine synthase, Beta ureidopropionase, BUP 1, BUP1, UPB1, ureidopropionase, beta
Observed MW:	43 kDa
Uniprot ID :	Q9UBR1

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:500-1:5000; IHC: 1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab09269(UPB1 Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab09269(UPB1 Antibody) at dilution of 1:1000