

anti- UQCR antibody

Product Information

Catalog No.:	FNab09278
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	ubiquinol-cytochrome c reductase, 6.4kDa subunit
Synonyms:	0710008D09Rik, Complex III subunit 10, Complex III subunit XI, Cytochrome b c1 complex subunit 10, QCR10, UQCR, UQCR11
Observed MW:	6 kDa
Uniprot ID :	O14957

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemical of paraffin-embedded human heart using FNab09278(UQCR antibody) at dilution of 1:50