

anti- UQCRC1 antibody

Product Information

Catalog No.:	FNab09280
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Immunogen information

Immunogen:	ubiquinol-cytochrome c reductase core protein I
Synonyms:	Complex III subunit 1, Core protein I, D3S3191, QCR1, UQCR1, UQCRC1
Observed MW:	53 kDa
UniprotID :	P31930

Application


Reactivity:	Human, Mouse, Rat, Monkey
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human heart using FNab09280(UQCRC1 antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab09280(UQCRC1 antibody) at dilution of 1:1000