

## anti- USE1 antibody

### Product Information

Catalog No.:	FNab09296
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER.

### Immunogen information

Immunogen:	unconventional SNARE in the ER 1 homolog(S. cerevisiae)
Synonyms:	MDS032, P31, SLT1, USE1, USE1 like protein, USE1L, Vesicle transport protein USE1
Observed MW:	29 kDa
Uniprot ID :	Q9NZ43

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:1000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09296(USE1 Antibody) at dilution of 1:300