

## anti- USF2 antibody

### Product Information

Catalog No.:	FNab09298
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Transcription factor that binds to a symmetrical DNA sequence(E-boxes)(5'-CACGTG-3') that is found in a variety of viral and cellular promoters.

### Immunogen information

Immunogen:	upstream transcription factor 2, c-fos interacting
Synonyms:	USF2, FOS interacting protein, bHLHb12, Upstream stimulatory factor 2, FIP
Observed MW:	44 kDa
Uniprot ID :	Q15853

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	

1


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09298(USF2 antibody) at dilution of 1:1000