

anti- VAPB antibody

Product Information

Catalog No.:	FNab09371
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9C9
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Participates in the endoplasmic reticulum unfolded protein response(UPR) by inducing ERN1/IRE1 activity. Involved in cellular calcium homeostasis regulation.

Immunogen information

Immunogen:	VAMP(vesicle-associated membrane protein)-associated protein B and C
Synonyms:	ALS8, VAMP associated protein B/C, VAMP B, VAMP B/VAMP C, VAMP C, VAP B, VAP B/VAP C, VAP C, VAPB
Observed MW:	27 kDa
Uniprot ID :	O95292

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab09371(VAPB Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab09371(VAPB Antibody) at dilution of 1:1000